



Parkinson's Application User- Manual

Version: 1.0

Last Updated: June 2025

Developed by: Sai Siddharth B S, Sourabh Mekkalike, Sripathi Baliga J, Vidya R

Introduction

This application is designed to help detect early signs of **Parkinson's Disease** through three non-invasive tests: **Speech Analysis**, **Handwriting Analysis**, and **Gait Analysis**. The app uses ML models trained on medical data to assist in early screening and monitoring.

System Requirements

- **Operating System:** Windows 10 or higher
- **RAM:** 8 GB minimum
- **Storage:** 1 GB free
- **Python:** Python 3.9 – 3.11 (only needed for script execution)
- **Optional Hardware:** Microphone, Writing pad, Webcam





Installation Instructions

1. Visit the official website or shared link.
2. Download the Parkinson_Detection.exe file.
3. Double-click to run (No installation required).
4. If prompted by Windows Defender, click **"More Info"** → **"Run"**
5. You're ready to begin!



How to Use the App


Main Menu Navigation

Choose from:



-  **Speech Analysis**
-  **Handwriting Analysis**
-  **Gait Analysis**
-  **Final Report Generator**

Speech Analysis

- Click **"Browse"** to select a .wav file of the patient reading a given paragraph.
- A waveform animation shows during analysis.
- After processing, results are shown as:
 -  **Healthy**
 -  **Parkinson Detected**

 *Tip: Use a quiet room and a good-quality microphone for best results.*

Handwriting Analysis

- Choose 4 spiral or pattern-drawing images using the provided radio buttons.
- Click "**Process**" to analyze them.
- Each image is individually evaluated and returns:
 -  **Healthy**
 -  **Parkinson Symptoms**

 *Tip: Scan clear, high-contrast images or use a drawing tablet if available.*

Gait Analysis

- Click "**Browse Gait Video**" and select a .mp4 or .avi file.
- The system analyzes steps, cadence, stride length, and posture angle.
- Output is displayed and used in the final report.

 *Tip: Use a well-lit video of the patient walking toward or away from the camera.*

Final Report

- Click "**Generate Report**" after completing the tests.
- Enter **Patient Name** and **ID**.

- The app generates a **customized diagnosis report** summarizing:
 - Speech results
 - Handwriting results
 - Gait analysis
 - Final assessment

Sample Output

=== PARKINSON'S DETECTION REPORT ===

Patient Name: John Doe

Patient ID: 1234

Date: 2025-05-23

HANDWRITING: 3/4 drawings indicate Parkinson

SPEECH: Parkinson Detected

GAIT: Abnormal cadence and stride length

FINAL ASSESSMENT:

 **STRONG INDICATORS OF PARKINSON'S DISEASE DETECTED**

Guidelines for Best Results

Test Type	Guideline
Speech	Use a quiet room, speak clearly into mic
Handwriting	Use scanned or tablet input with visible lines
Gait	Record video from the side/front, with good light

Privacy and Security

- All data is processed **locally** on your system
- No data is uploaded anywhere
- **100% offline** and private

Troubleshooting

Problem	Solution
App doesn't open	Right-click > Run as Administrator
.wav file not accepted	Ensure it's in .wav format only
Gait video too short	Use at least a 10-second walking video

Future Features

- Patient history tracking
- Medication response logging
- Online doctor dashboard integration
- Progress graphing
- Mobile app release